



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Drop-n-Stain EverBrite(TM) Mounting Medium with DAPI BOT-23009-T

Article Name	Drop-n-Stain EverBrite(TM) Mounting Medium with DAPI
Biozol Catalog Number	BOT-23009-T
Supplier Catalog Number	23009-T
Alternative Catalog Number	BOT-23009-T
Manufacturer	Biotium
Category	Zellen/Zellkultur
Product Description	Drop-n-Stain EverBrite(TM) Mounting Medium with DAPI is wet-set antifade in a convenient dropper bottle for one step mounting and counterstaining. EverBrite(TM) protects fluorescent dyes from photobleaching, and is compatible with cyanine-based dyes
Source	Synthetic