



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Green Downward cADDis cAMP Assay - 60mL MTM-D0200G-060

Article Name	Green Downward cADDis cAMP Assay - 60mL
Biozol Catalog Number	MTM-D0200G-060
Supplier Catalog Number	D0200G-060
Alternative Catalog Number	MTM-D0200G-060
Manufacturer	Montana Molecular
Category	Kits/Assays
Product Description	Measure cyclic AMP (cAMP) signaling kinetics, screen GPCRs for Gs activity, and more with the live cell, fluorescent, cADDis cAMP assay. Add the cADDis sensor to your cells and measure robust signals the next day on an automated plate reader or fluor