



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Green Downward DAG Assay - 30mL MTM-D0300G-030

Article Name	Green Downward DAG Assay - 30mL
Biozol Catalog Number	MTM-D0300G-030
Supplier Catalog Number	D0300G-030
Alternative Catalog Number	MTM-D0300G-030
Manufacturer	Montana Molecular
Category	Kits/Assays
Product Description	Measure diacylglycerol (DAG) signaling kinetics, screen GPCRs for Gq activity, and more with this live cell fluorescent assay. Add the DAG sensor to your cells and measure robust signals the next day on an automated plate reader or fluorescence imagi