

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Borealis Arrestin V2R Assay - 30mL** **MTM-D2000G-030**

|                            |   |
|----------------------------|---|
| Article Name               | Borealis Arrestin V2R Assay - 30mL  |
| Biozol Catalog Number      | MTM-D2000G-030  |
| Supplier Catalog Number    | D2000G-030  |
| Alternative Catalog Number | MTM-D2000G-030  |
| Manufacturer               | Montana Molecular   |
| Category                   | Kits/Assays   |
| Product Description        | Measure arrestin signaling kinetics in real time with the live cell, fluorescent Borealis arrestin - Vasopressin Receptor (V2R) assay. Add the unmodified V2R receptor and Borealis sensor to cells and measure robust signals the next day on an automate... |