



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Red Nucleus-targeted cADDis cAMP Assay - 30mL MTM-U0221R-030

Article Name	Red Nucleus-targeted cADDis cAMP Assay - 30mL
Biozol Catalog Number	MTM-U0221R-030
Supplier Catalog Number	U0221R-030
Alternative Catalog Number	MTM-U0221R-030
Manufacturer	Montana Molecular
Category	Kits/Assays
Product Description	Measure real-time cAMP signaling kinetics in the nucleus cilia with this live cell fluorescent assay. Add NLS-tagged cADDis to your cells and measure robust signals the next day on a fluorescence imaging system. Increasing red fluorescence indicates