



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Green Gi cADDis cAMP Assay - 30mL MTM-X0200G-030

Article Name	Green Gi cADDis cAMP Assay - 30mL
Biozol Catalog Number	MTM-X0200G-030
Supplier Catalog Number	X0200G-030
Alternative Catalog Number	MTM-X0200G-030
Manufacturer	Montana Molecular
Category	Kits/Assays
Product Description	Measure cyclic AMP (cAMP) signaling kinetics, screen GPCRs for Gi activity, and more with the live cell, fluorescent, cADDis cAMP assay. Add the cADDis sensor to your cells and measure robust signals the next day on an automated plate reader or fluor