

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Green Gi cADDIs cAMP Assay - 60mL **MTM-X0200G-060**

| | |
|----------------------------|---|
| Article Name | Green Gi cADDIs cAMP Assay - 60mL |
| Biozol Catalog Number | MTM-X0200G-060 |
| Supplier Catalog Number | X0200G-060 |
| Alternative Catalog Number | MTM-X0200G-060 |
| Manufacturer | Montana Molecular |
| Category | Kits/Assays |
| Product Description | Measure cyclic AMP (cAMP) signaling kinetics, screen GPCRs for Gi activity, and more with the live cell, fluorescent, cADDIs cAMP assay. Add the cADDIs sensor to your cells and measure robust signals the next day on an automated plate reader or fluor... |