

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Recombinant Enterobacteria phage T4 Single-stranded DNA-binding protein (32) CSB-YP022706EDZ

Recombinant Enterobacteria phage T4 Single-stranded DNA-binding protein (32)
CSB-YP022706EDZ
CSB-YP022706EDZ
CSB-YP022706EDZ-1, CSB-YP022706EDZ-100, CSB-YP022706EDZ-20
Cusabio
Proteine/Peptide
35.0 kDa
C-terminal 6xHis-tagged
P03695
If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Yeast
1-301aa
Greater than 85% as determined by SDS-PAGE.
Liquid or Lyophilized powder

Sequence

MFKRKSTAELAAQMAKLNGNKGFSSEDKGEWKLKLDNAGNGQAVIRFLPSK NDEQAPFAILVNHGFKKNGKWYIETCSSTHGDYDSCPVCQYISKNDLYNTDN KEYSLVKRKTSYWANILVVKDPAAPENEGKVFKYRFGKKIWDKINAMIAVDVE MGETPVDVTCPWEGANFVLKVKQVSGFSNYDESKFLNQSAIPNIDDESFQKE LFEQMVDLSEMTSKDKFKSFEELNTKFGQVMGTAVMGGAAATAAKK