

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

ELISA Kit for Nicotinic Acid Adenine Dinucleotide Phosphate (NAADP) USC-CEG366GE

Article Name	ELISA Kit for Nicotinic Acid Adenine Dinucleotide Phosphate (NAADP)
Biozol Catalog Number	USC-CEG366GE
Supplier Catalog Number	CEG366Ge
Alternative Catalog Number	USC-CEG366GE-48, USC-CEG366GE-96
Manufacturer	Cloud-Clone
Category	Kits/Assays
Species Reactivity	All
Range	123.5-10,000ng/mL
Sensitivity	The minimum detectable dose of this kit is typically less than 48.7ng/mL
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Application Notes	Enzyme-linked immunosorbent assay for Antigen Detection.