

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

GFP/CHO-K1 Stable Cell Line ABI-14-901ACL

Article Name	GFP/CHO-K1 Stable Cell Line
Biozol Catalog Number	ABI-14-901ACL
Supplier Catalog Number	14-901ACL
Alternative Catalog Number	ABI-14-901ACL-1VIAL
Manufacturer	Abeomics
Category	Zellen/Zellkultur
Application	FA, FACS
Product Description	GFP/CHO-K1 Stable Cell Line is a stably transfected CHO-K1 cell line which expresses enhanced green fluorescent protein (eGFP). Sequence data: Amino acid sequence of eGFP MVSKEELFTGVVPIVLDGVDVNGHKFSVSGEGEGDATYGLTLKFIC TTGKLPVPWPTLVTTLTLYGVQCFSRYPDH...
Application Dilute	Application:.. Screen for GFP through Flow Cytometry. Screen for GFP through Fluorescence Microscopy. Culture conditions: Cells should be grown at 37°C with 5% CO ₂ using DMEM medium (w/ L-Glutamine, 4.5g/L Glucose and Sodium Pyruvate) supplemented with 10