

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Monosan Plus (AP) Anti-Mouse (Fast Red)Kit, 1 Kit (8 ml/80 tests) MON-APP119

Article Name	Monosan Plus (AP) Anti-Mouse (Fast Red)Kit, 1 Kit (8 ml/80 tests)
Biozol Catalog Number	MON-APP119
Supplier Catalog Number	MON-APP119
Alternative Catalog Number	MON-APP119
Manufacturer	Monosan
Category	Kits/Assays
Application	IF, IHC-Fr, IHC-P
Product Description	The purpose of immunohistochemical staining is to make tissue and cell antigens visible. The Plus AP Kit, Mouse is a highly sensitive detection kit intended for use in immunohistochemistry and immunocytochemistry. The method is based on the streptavi