



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Monosan Wash Buffer (20x), 2,500 ml (for 50 liter) MON-APP166

Article Name	Monosan Wash Buffer (20x), 2,500 ml (for 50 liter)
Biozol Catalog Number	MON-APP166
Supplier Catalog Number	MON-APP166
Alternative Catalog Number	MON-APP166
Manufacturer	Monosan
Category	Sonstiges
Application	IF, IHC-Fr, IHC-P
Product Description	Immunohistochemical staining procedures consist of sequential incubation steps with blocking solutions, antibodies and secondary reagents, enzymes and chromogenic substrates carried out on tissue sections. Washing away the applied reagents after each