

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Monosan Fast (AP) One-Step Polymer anti-Mouse/Rabbit/Rat, 6 ml/60 tests MON-APP207

Article Name	Monosan Fast (AP) One-Step Polymer anti-Mouse/Rabbit/Rat, 6 ml/60 tests
Biozol Catalog Number	MON-APP207
Supplier Catalog Number	MON-APP207
Alternative Catalog Number	MON-APP207
Manufacturer	Monosan
Category	Sonstiges
Application	IF, IHC-Fr, IHC-P
Product Description	The purpose of the immunohistochemical staining is to make tissue and cell antigens visible. The Fast AP One-Step Polymer anti-Mouse/Rabbit is a highly sensitive detection reagent intended for use in immunohistochemistry and immunocytochemistry. The