

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

BrightVision, 2 steps detection system goat anti-mouse/rabbit AP MON-APP904

| Article Name | BrightVision, 2 steps detection system goat anti-mouse/rabbit AP |
|----------------------------|--|
| Biozol Catalog Number | MON-APP904 |
| Supplier Catalog Number | MON-APP904 |
| Alternative Catalog Number | MON-APP904 |
| Manufacturer | Monosan |
| Category | Sonstiges |
| Application | IHC-P |
| Product Description | The BrightVision detection system peroxidase Goat Anti-Mouse/Rabbit AP, is a Ready-to-Use system that has been manufactured to give an optimal staining, when using the protocol advised in this IFU. Prior to staining some routine fixed, paraffinembe |