

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## BrightVision, colour 2 steps detection system goat anti-mouse/rabbit AP MON-APP905

Article Name	BrightVision, colour 2 steps detection system goat anti-mouse/rabbit AP
Biozol Catalog Number	MON-APP905
Supplier Catalog Number	MON-APP905
Alternative Catalog Number	MON-APP905
Manufacturer	Monosan
Category	Sonstiges
Application	IHC-P
Product Description	The BrightVision detection system peroxidase Goat Anti- Mouse/Rabbit AP, is a Ready-to-Use system that has been manufactured to give an optimal staining, when using the protocol advised in this IFU. Prior to staining some routine fixed, paraffin- embe