

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

MGMT (Marker of Suppressor of Mutation and Carcinogenesis) Antibody, Clone: [MGMT/8186R], Rabbit, Monoclonal NBT-4255-RBM8-P1

Article Name	MGMT (Marker of Suppressor of Mutation and Carcinogenesis) Antibody, Clone: [MGMT/8186R], Rabbit, Monoclonal
Biozol Catalog Number	NBT-4255-RBM8-P1
Supplier Catalog Number	4255-RBM8-P1
Alternative Catalog Number	NBT-4255-RBM8-P1-100
Manufacturer	NeoBiotechnologies
Host	Rabbit
Category	Antikörper
Application	IHC
Species Reactivity	Human
Immunogen	Recombinant fragment (around aa1-100) of human MGMT (exact sequence is proprietary)
Product Description	Cancer chemotherapeutic alkylating agents (e.g. BCNU,) act by inducing formation of lethal cross links at the O6 alkylguanine position in DNA. MGMT transfers alkyl adducts from the O6 position of guanine in DNA (prior to cross link formation) to a c
Clonality	Monoclonal
Clone Designation	[MGMT/8186R]
Molecular Weight	21.646 kDa
NCBI	4255

UniProt	P16455
Form	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double Concentrate by Protein A/G. Protein