

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anti-Human CD3-FITC/CD19-APC/CD16+CD56-PE Cocktail, Clone: [OKT-3CB193G85.1H11], Monoclonal ELA-E-AB-FC0007

Article Name	Anti-Human CD3-FITC/CD19-APC/CD16+CD56-PE Cocktail, Clone: [OKT-3CB193G85.1H11], Monoclonal
Biozol Catalog Number	ELA-E-AB-FC0007
Supplier Catalog Number	E-AB-FC0007
Alternative Catalog Number	ELA-E-AB-FC0007-100, ELA-E-AB-FC0007-20, ELA-E-AB-FC0007-100X2
Manufacturer	Elabscience
Category	Antikörper
Application	FC
Species Reactivity	Human
Clonality	Monoclonal
Concentration	5 uL/Test
Clone Designation	[OKT-3CB193G85.1H11]
Buffer	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Application Notes	For whole blood samples, add 5 uL Anti-Human CD3-FITC/CD19-APC/CD16+CD56-PE Cocktail to 100 uL anticoagulant-treated blood sample. Mix and incubate the sample at 4C in the dark for 30 min. Remove red blood cells with RBC lysis solution following the man