

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Anti-Human CD3-FITC/CD4-PE/CD8a-PE/Cyanine7/CD45-PerCP Cocktail, Clone: [OKT-3RPA-T4OKT-8HI30], Monoclonal ELA-E-AB-FC0009**

Article Name	Anti-Human CD3-FITC/CD4-PE/CD8a-PE/Cyanine7/CD45-PerCP Cocktail, Clone: [OKT-3RPA-T4OKT-8HI30], Monoclonal
Biozol Catalog Number	ELA-E-AB-FC0009
Supplier Catalog Number	E-AB-FC0009
Alternative Catalog Number	ELA-E-AB-FC0009-100, ELA-E-AB-FC0009-20, ELA-E-AB-FC0009-100X2
Manufacturer	Elabscience
Category	Antikörper
Application	FC
Species Reactivity	Human
Clonality	Monoclonal
Concentration	5 uL/Test
Clone Designation	[OKT-3RPA-T4OKT-8HI30]
Buffer	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Application Notes	For whole blood samples, add 5 uL Anti-Human CD3-FITC/CD4-PE/CD8a-PE/Cyanine7/CD45-PerCP Cocktail to 100 uL anticoagulant-treated blood sample. Mix and incubate the sample at 4C in the dark for 30 min. Remove red blood cells with RBC lysis solution foll