

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Anti-Human CD8a-FITC/CD4-PE/CD3-PE/Cyanine7/CD45-PerCP Cocktail, Clone: [OKT-8RPA-T4OKT-3HI30], Monoclonal ELA-E-AB-FC0010

Article Name	Anti-Human CD8a-FITC/CD4-PE/CD3-PE/Cyanine7/CD45-PerCP Cocktail, Clone: [OKT-8RPA-T40KT-3HI30], Monoclonal
Biozol Catalog Number	ELA-E-AB-FC0010
Supplier Catalog Number	E-AB-FC0010
Alternative Catalog Number	ELA-E-AB-FC0010-100, ELA-E-AB-FC0010-20, ELA-E-AB-FC0010- 100X2
Manufacturer	Elabscience
Category	Antikörper
Application	FC
Species Reactivity	Human
Clonality	Monoclonal
Concentration	5 uL/Test
Clone Designation	[OKT-8RPA-T4OKT-3HI30]
Buffer	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Application Notes	For whole blood samples, add 5 uL Anti-Human CD8a-FITC/CD4-PE/CD3-PE/Cyanine7/CD45-PerCP Cocktail to 100 uL anticoagulant-treated blood sample. Mix and incubate the sample at 4C in the dark for 30 min. Remove red blood cells with RBC lysis solution foll