

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

SARS-CoV-2 Variants Total IgG 6-plex Panel (Flow Cytometry Multiplex Bead Assay, Non-premixed) **ABS-FCM-B08R-96TESTS**

Article Name	SARS-CoV-2 Variants Total IgG 6-plex Panel (Flow Cytometry Multiplex Bead Assay, Non-premixed)
Biozol Catalog Number	ABS-FCM-B08R-96TESTS
Supplier Catalog Number	FCM-B08R-96tests
Alternative Catalog Number	ABS-FCM-B08R-96TESTS
Manufacturer	AcroBiosystems
Category	Kits/Assays
Species Reactivity	Virus
Product Description	Since December 2019, the emergence of severe acute respiratory syndrome coronavirus 2 (SARSCoV2) and its associated disease, COVID19, has caused a devastating pandemic worldwide.As the virus spreads globally, the continuous emergence of new mutant st...
Source	Since December 2019, the emergence of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and its associated disease, COVID-19, has caused a devastating pandemic worldwide.As the virus spreads globally, the continuous emergence of new mutant str
Form	Solid
Target	Spike RBD
Application Notes	The kit is developed for qualitative detection of Anti-SARS-CoV-2 Variants (Alpha, Beta, Gamma, Delta, Omicron) total IgG antibody in human serum.