

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### AP-TNAP/ALPL Polyclonal Antibody, PerCP Conjugated, Rabbit BSS-BS-1535R-PERCP

Article Name	AP-TNAP/ALPL Polyclonal Antibody, PerCP Conjugated, Rabbit
Biozol Catalog Number	BSS-BS-1535R-PERCP
Supplier Catalog Number	BS-1535R-PERCP
Alternative Catalog Number	BSS-BS-1535R-PERCP-100,BSS-BS-1535R-PERCP-20
Manufacturer	Bioss
Host	Rabbit
Category	Antikörper
Application	IF, WB
Species Reactivity	Human, Mouse
Immunogen	52-150/524
Conjugation	PerCP
Product Description	Alkaline phosphatase (ALP) removes phosphate groups from the 5 end of DNA and RNA, and from proteins, at high pH. Most mammals have 4 different isozymes: placental, placental like, intestinal and non tissue specific (found in liver, kidney and bone)....
Clonality	Polyclonal
Concentration	1ug/ul
NCBI	<a href="#">11647</a>
Buffer	Aqueous buffered solution containing 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Source	KLH conjugated synthetic peptide derived from AP-TNAP/ALPL
Purity	Purified by Protein A.
Target	AP-TNAP/ALPL
Application Dilute	WB(1:300-5000), IF(IHC-P)(1:50-200), IF(IHC-F)(1:50-200), IF(ICC)(1:50-200)