



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## MAPRE3 Peptide PRS-42-280P

Article Name	MAPRE3 Peptide
Biozol Catalog Number	PRS-42-280P
Supplier Catalog Number	42-280P
Alternative Catalog Number	PRS-42-280P-0.1
Manufacturer	ProSci
Category	Proteine/Peptide
Application	Blocking
NCBI	22924
UniProt	Q9UPY8
Buffer	MAPRE3 peptide is supplied as 100 $\mu g$ of dried peptide. When the peptide is reconstituted in 200 $\mu L$ water, the concentration will be 0.5 mg/mL.
Form	Lyophilized
Application Notes	MAPRE3 Peptide is used for blocking the activity of MAPRE3 antibody (catalog number 42-280). The MAPRE3 antibody ships at 0.5 mg/mL and the MAPRE3 Peptide ships as a 100 ug pellet. When peptide is reconstituted in 200 uL water, the concentration will be