



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Sap30L Peptide PRS-42-835P

Article Name	Sap30L Peptide
Biozol Catalog Number	PRS-42-835P
Supplier Catalog Number	42-835P
Alternative Catalog Number	PRS-42-835P-0.1
Manufacturer	ProSci
Category	Proteine/Peptide
Application	Blocking
NCBI	50724
UniProt	Q5SQF8
Buffer	Sap30L peptide is supplied as 100 $\mu g$ of dried peptide. When the peptide is reconstituted in 200 $\mu L$ water, the concentration will be 0.5 mg/mL.
Form	Lyophilized
Application Notes	Sap30L Peptide is used for blocking the activity of Sap30L antibody (catalog number 42-835). The Sap30L antibody ships at 0.5 mg/mL and the Sap30L Peptide ships as a 100 $\mu$ g pellet. When peptide is reconstituted in 200 $\mu$ L water, the concentration will be