



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## TIMP1 Peptide PRS-42-914P

| Article Name               | TIMP1 Deatide   |
|----------------------------|---|
| Article Name               | TIMP1 Peptide   |
| Biozol Catalog Number      | PRS-42-914P   |
| Supplier Catalog Number    | 42-914P   |
| Alternative Catalog Number | PRS-42-914P-0.1   |
| Manufacturer               | ProSci  |
| Category                   | Proteine/Peptide  |
| Application                | Blocking  |
| NCBI                       | 7076  |
| UniProt                    | P01033  |
| Buffer                     | TIMP1 peptide is supplied as 100 $\mu g$ of dried peptide. When the peptide is reconstituted in 200 $\mu L$ water, the concentration will be 0.5 mg/mL.   |
| Form                       | Lyophilized   |
| Application Notes          | TIMP1 Peptide is used for blocking the activity of TIMP1 antibody (catalog number 42-914). The TIMP1 antibody ships at 0.5 mg/mL and the TIMP1 Peptide ships as a 100 $\mu$ g pellet. When peptide is reconstituted in 200 $\mu$ L water, the concentration will be 0.5 |