



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## KCNIP3 Peptide PRS-43-040P

| Article Name               | KCNIP3 Peptide   |
|----------------------------|--|
| Biozol Catalog Number      | PRS-43-040P  |
| Supplier Catalog Number    | 43-040P  |
| Alternative Catalog Number | PRS-43-040P-0.1  |
| Manufacturer               | ProSci   |
| Category                   | Proteine/Peptide   |
| Application                | Blocking   |
| NCBI                       | 30818  |
| UniProt                    | Q9Y2W7   |
| Buffer                     | KCNIP3 Peptide is supplied as 100 $\mu g$ of dried peptide. When the peptide is reconstituted in 200 $\mu L$ water, the concentration will be 0.5 mg/mL.   |
| Form                       | Lyophilized  |
| Application Notes          | KCNIP3 Antibody ships at 0.5 mg/mL and the KCNIP3 Peptide ships as a 100 $\mu$ g pellet. When peptide is reconstituted in 200 $\mu$ L water, the concentration will be 0.5 mg/mL, identical to that of the antibody. The best starting ratio would be 1:1 (which means m |