



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

perilipin 3 / TIP47 Peptide PRS-43-202P

Article Name	perilipin 3 / TIP47 Peptide
Biozol Catalog Number	PRS-43-202P
Supplier Catalog Number	43-202P
Alternative Catalog Number	PRS-43-202P-0.1
Manufacturer	ProSci
Category	Proteine/Peptide
Application	Blocking
NCBI	10226
UniProt	O60664
Buffer	perilipin 3 / TIP47 peptide is supplied as 100 μg of dried peptide. When the peptide is reconstituted in 200 μL water, the concentration will be 0.5 mg/mL.
Form	Lyophilized
Application Notes	perilipin 3 / TIP47 peptide is used for blocking the activity of perilipin 3 / TIP47 antibody (catalog number 43-202). The perilipin 3 / TIP47 antibody ships at 0.5 mg/mL and the STK36 Peptide ships as a 100 ug pellet. When peptide is reconstituted in 2