

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

DACH2 Peptide PRS-46-650P

| | |
|----------------------------|---|
| Article Name | DACH2 Peptide |
| Biozol Catalog Number | PRS-46-650P |
| Supplier Catalog Number | 46-650P |
| Alternative Catalog Number | PRS-46-650P-0.1 |
| Manufacturer | ProSci |
| Category | Proteine/Peptide |
| Application | Blocking |
| NCBI | 117154 |
| UniProt | Q96NX9 |
| Buffer | DACH2 peptide is supplied as 100 µg of dried peptide. When the peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL. |
| Form | Lyophilized |
| Application Notes | DACH2 Peptide is used for blocking the activity of DACH2 antibody (catalog number 46-650). The DACH2 antibody ships at 0.5 mg/mL and the DACH2 Peptide ships as a 100 ug pellet. When peptide is reconstituted in 200 uL water, the concentration will be 0.5 |