



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Frizzled 4 Peptide PRS-46-660P

Article Name	Frizzled 4 Peptide
Biozol Catalog Number	PRS-46-660P
Supplier Catalog Number	46-660P
Alternative Catalog Number	PRS-46-660P-0.1
Manufacturer	ProSci
Category	Proteine/Peptide
Application	Blocking
NCBI	8322
UniProt	Q9ULV1
Buffer	Frizzled 4 peptide is supplied as 100 $\mu g$ of dried peptide. When the peptide is reconstituted in 200 $\mu L$ water, the concentration will be 0.5 mg/mL.
Form	Lyophilized
Application Notes	Frizzled 4 Peptide is used for blocking the activity of Frizzled 4 antibody (catalog number 46-660). The Frizzled 4 antibody ships at 0.5 mg/mL and the Frizzled 4 Peptide ships as a 100 ug pellet. When peptide is reconstituted in 200 uL water, the concen