



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

SAPAP Detection Set, Polyclonal PRS-PSI-1826

Article Name	SAPAP Detection Set, Polyclonal
Biozol Catalog Number	PRS-PSI-1826
Supplier Catalog Number	PSI-1826
Alternative Catalog Number	PRS-PSI-1826-1
Manufacturer	ProSci
Category	Antikörper
Application	IF, IHC, WB
Species Reactivity	Human
Immunogen	Rabbit polyclonal antibodies were raised against peptides corresponding to amino acid sequences from each of the corresponding proteins.
Clonality	Polyclonal
Concentration	Antibody 1 mg/mL
Buffer	PBS containing 0.02% sodium azide.
Form	Liquid
Application Dilute	Optimal dilutions for each application to be determined by the researcher.
Application Notes	These polyclonal antibodies can be used for detection of SAPAP1 - 4 by immunoblot at 0.5 - 2 μ g/mL, and for detection of SAPAP1 - 3 by immunohistochemistry at 2.5 - 5 μ g/mL, and Immunoflourescence.