

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Loading Control Detection Set Alpha-Tubulin, Beta-Actin, GAPDH, Polyclonal PRS-PSI-1836

Article Name	Loading Control Detection Set Alpha-Tubulin, Beta-Actin, GAPDH, Polyclonal
Biozol Catalog Number	PRS-PSI-1836
Supplier Catalog Number	PSI-1836
Alternative Catalog Number	PRS-PSI-1836-0.02, PRS-PSI-1836-0.1
Manufacturer	ProSci
Category	Antikörper
Application	WB
Species Reactivity	Human
Immunogen	Rabbit polyclonal antibodies were raised against peptides corresponding to amino acid sequences from each of the corresponding proteins.
Clonality	Polyclonal
Concentration	Antibody 1 mg/mL
Buffer	PBS containing 0.02% sodium azide.
Form	Liquid
Application Dilute	Optimal dilutions for each application to be determined by the researcher.