

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

PD-1 Mapping Peptide Set PRS-RF16000-02

| | |
|----------------------------|--|
| Article Name | PD-1 Mapping Peptide Set |
| Biozol Catalog Number | PRS-RF16000-02 |
| Supplier Catalog Number | RF16000-02 |
| Alternative Catalog Number | PRS-RF16000-02-0.05 |
| Manufacturer | ProSci |
| Category | Proteine/Peptide |
| Application | Blocking |
| Concentration | 200 ug/mL |
| NCBI | 5133 |
| UniProt | Q15116 |
| Buffer | PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide |
| Form | Liquid |
| Application Notes | PD-1 mapping peptides can be used to map the epitopes of PD-1 antibodies. |