

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## TWEAK Antibody, Unconjugated, Goat, Polyclonal PRS-XP-5290

Article Name	TWEAK Antibody, Unconjugated, Goat, Polyclonal
Biozol Catalog Number	PRS-XP-5290
Supplier Catalog Number	XP-5290
Alternative Catalog Number	PRS-XP-5290-0.1
Manufacturer	ProSci
Host	Goat
Category	Antikörper
Application	ELISA, WB
Species Reactivity	Human
Immunogen	Produced from sera of goats pre-immunized with highly pure (>98%) recombinant hTWEAK (human TWEAK).
Conjugation	Unconjugated
Clonality	Polyclonal
Concentration	batch dependent
NCBI	8742
UniProt	O43508
Form	Lyophilized
Application Dilute	Centrifuge vial prior to opening.

**Application Notes** 

ELISA:To detect hTWEAK by direct ELISA (using 100  $\mu$ L/well antibody solution) a concentration of at least 0.5  $\mu$ g/mL of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the de