

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## TWEAK Antibody (biotin), Biotin, Goat, Polyclonal PRS-XP-5290BT

Article Name	TWEAK Antibody (biotin), Biotin, Goat, Polyclonal
Biozol Catalog Number	PRS-XP-5290BT
Supplier Catalog Number	XP-5290Bt
Alternative Catalog Number	PRS-XP-5290BT-0.05
Manufacturer	ProSci
Host	Goat
Category	Antikörper
Application	ELISA, WB
Species Reactivity	Human
Immunogen	Produced from sera of goats pre-immunized with highly pure (>98%) recombinant hTWEAK (human TWEAK).
Conjugation	Biotin
Clonality	Polyclonal
Concentration	batch dependent
NCBI	8742
UniProt	O43508
Form	Lyophilized
Application Dilute	Centrifuge vial prior to opening.

**Application Notes** 

ELISA:Sandwich:To detect hTWEAK by sandwich ELISA (using 100  $\mu\text{L/well}$  antibody solution) a concentration of 0.25 - 1.0  $\mu\text{g/mL}$  of this antibody is required. This biotinylated polyclonal antibody, in conjunction with our Polyclonal Anti-Human TWEAK (