

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Porin P Protein, Pseudomonas aeruginosa, Recombinant (His & SUMO) TGM-TMPH-03175

Article Name	Porin P Protein, Pseudomonas aeruginosa, Recombinant (His & SUMO)
Biozol Catalog Number	TGM-TMPH-03175
Supplier Catalog Number	TMPH-03175
Alternative Catalog Number	TGM-TMPH-03175-1MG,TGM-TMPH-03175-20UG,TGM-TMPH-03175- 100UG
Manufacturer	TargetMol
Category	Biochemikalien
Product Description	Anion specific, the binding site has higher affinity for phosphate than chloride ions. Porin O has a higher affinity for polyphosphates (tripolyphosphate and pyrophosphate) while porin P has a higher affinity for orthophosphate. Porin P Protein, Pseu
Molecular Weight	61.2 kDa (predicted)