

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mycoplasma PCR Detection Kit ASC-BMC1040

| | |
|----------------------------|---|
| Article Name | Mycoplasma PCR Detection Kit |
| Biozol Catalog Number | ASC-BMC1040 |
| Supplier Catalog Number | BMC1040 |
| Alternative Catalog Number | ASC-BMC1040-200, ASC-BMC1040-1000 |
| Manufacturer | Abbkine Scientific |
| Category | Kits/Assays |
| Product Description | This Mycoplasma PCR Detection Kit utilizes Nested PCR technology to amplify and detect the conserved region specific fragments of Mycoplasma 16S-23S rRNA sequence. It can be directly used for rapid detection of biomaterial samples such as cell cultur... |
| Tag | Mycoplasma |
| Form | Liquid solution |
| Samples | Cell culture supernatant, laboratory animal secretions, serum and other biomaterial samp |
| Application Notes | In order to prevent contamination of the Positive Control Template, it is necessary to add the Positive Control Template to the positive control group after adding samples in the test group and the negative control group. |