

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red) ASC-BMU109-EN

Article Name	SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red)
Biozol Catalog Number	ASC-BMU109-EN
Supplier Catalog Number	BMU109-EN
Alternative Catalog Number	ASC-BMU109-EN-100
Manufacturer	Abbkine Scientific
Category	Sonstiges
Application	Cell Culture
Product Description	SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red) contains 0.25% trypsin and 0.06% EDTA, and the pH value is 7.2-7.4. After filtration and sterilization, the digestive solution can be directly used for the digestion of cultured cells or so...
Tag	Trypsin-EDTA Solution
Form	Liquid solution
Samples	Cells and Tissues
Application Notes	In the process of using SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red), special attention should be paid to avoid bacterial contamination of the digest. The digestion time of trypsin cell digestive solution should not be too long, otherw