

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **TAP-Tag (Mix-mA?) Monoclonal Antibody, AbBy Fluor 405 Conjugated, Clone: [Mix-mA?], AbFluor 405, Mouse BSS-BSM-33133M-BF405**

Article Name	TAP-Tag (Mix-mA?) Monoclonal Antibody, AbBy Fluor 405 Conjugated, Clone: [Mix-mA?], AbFluor 405, Mouse
Biozol Catalog Number	BSS-BSM-33133M-BF405
Supplier Catalog Number	BSM-33133M-BF405
Alternative Catalog Number	BSS-BSM-33133M-BF405-100
Manufacturer	Bioss
Host	Mouse
Category	Antikörper
Application	WB
Conjugation	AbFluor 405
Product Description	Tandem affinity purification (TAP) is a purification technique for studying proteinprotein interactions. It involves creating a fusion protein with a designed piece, the TAP tag, on the end. The original TAP method involves the fusion of the TAP tag ...
Clonality	Monoclonal
Concentration	1ug/ul
Clone Designation	[Mix-mA?]
Buffer	Aqueous buffered solution containing 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Source	Recombinant TAP-Tag

Purity	Purified by Protein A.
Target	TAP-Tag
Application Dilute	WB(1:300-5000)