



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Proteinase K (25X) GTX00761

Article Name	Proteinase K (25X)
Biozol Catalog Number	GTX00761
Supplier Catalog Number	GTX00761
Alternative Catalog Number	GTX00761-2
Manufacturer	GeneTex
Category	Sonstiges
Application	IHC-P
Product Description	Proteinase K is specifically formulated to achieve proteolytic digestion of formalin fixed, paraffin embedded tissue sections prior to performing in situ hybridization procedures. Proteolytic digestion of formalin fixed tissues improves the accessibi
Form	Liquid
Application Notes	Preparation of working solutions : 1. Place 100µl Proteinase K solution in a glass tube, and 2.4ml distilled water and mix. Protocol Recommendations : 1. Optimal digestion conditions will vary depending on fixation conditions and should be op