

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Trident Blue Prestained Protein Ladder GTX16376

Article Name	Trident Blue Prestained Protein Ladder
Biozol Catalog Number	GTX16376
Supplier Catalog Number	GTX16376
Alternative Catalog Number	GTX16376-500
Manufacturer	GeneTex
Category	Sonstiges
Application	WB
Form	Liquid
Application Notes	3 µl or 5 µl per loading for clear visualization during electrophoresis on 15-well or 10-well mini-gel, respectively. 1.5~2.5 µl per well for general Western transferring. Apply more for thicker (> 1.5 mm) or larger gel.