

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Nestin antibody, Unconjugated, Gallus, Polyclonal GTX85457

Article Name	Nestin antibody, Unconjugated, Gallus, Polyclonal
Biozol Catalog Number	GTX85457
Supplier Catalog Number	GTX85457
Alternative Catalog Number	GTX85457-50
Manufacturer	GeneTex
Host	Gallus
Category	Antikörper
Application	ICC, IHC, WB
Species Reactivity	Mouse
Immunogen	Peptides
Conjugation	Unconjugated
Clonality	Polyclonal
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Molecular Weight	177
Sensitivity	Nestin
NCBI	10763
UniProt	P48681
Buffer	PBS, 0.02% Sodium azide.

Form	Liquid
Application Notes	Recommended starting dilutions: - For IHC: Use at 1:3000 - 1:6000. - For ICC/IF:Use at 1:3000 - 1:6000. - For WB: Use at 1:5000 - 1:10000. Not yet tested in other applications. Optimal working dilutions should be determined experimentally by th