

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

SUMO2 + SUMO3 (cleaved Gly93) antibody, Unconjugated, Rabbit, Polyclonal GTX86908

Article Name	SUMO2 + SUMO3 (cleaved Gly93) antibody, Unconjugated, Rabbit, Polyclonal
Biozol Catalog Number	GTX86908
Supplier Catalog Number	GTX86908
Alternative Catalog Number	GTX86908-100
Manufacturer	GeneTex
Host	Rabbit
Category	Antikörper
Application	WB
Species Reactivity	Human
Immunogen	The antiserum was produced against synthesized peptide derived from human SUMO2 + SUMO3 (44-93).
Conjugation	Unconjugated
Clonality	Polyclonal
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Sensitivity	SUMO2 + SUMO3 (cleaved Gly93) antibody detects endogenous levels of fragment of activated SUMO2/3 resulting from cleavage adjacent to Gly93.
NCBI	6612

UniProt	P61956, P55854
Buffer	PBS (without Mg and Ca) pH7.4, 150mM NaCl, 50% Glycerol, 0.02% Sodium azide.
Form	Liquid
Application Notes	WB: 1:500~1:1000. *Optimal dilutions/concentrations should be determined by the researcher.Not tested in other applications.