



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

PDIA2 / PDIP blocking peptide GTX89110-PEP

Article Name	PDIA2 / PDIP blocking peptide
Biozol Catalog Number	GTX89110-PEP
Supplier Catalog Number	GTX89110-PEP
Alternative Catalog Number	GTX89110-PEP-100
Manufacturer	GeneTex
Category	Proteine/Peptide
Species Reactivity	Human
Sensitivity	The peptide was used in the production of GTX89110.
NCBI	64714
Buffer	Reconstitute with 200 μ l distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H2O, no preservative.
Form	Lyophilized powder
Sequence	C-NLETTKKYAPVD