



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## PI3 kinase p110 alpha blocking peptide GTX89196-PEP

Article Name	PI3 kinase p110 alpha blocking peptide
Biozol Catalog Number	GTX89196-PEP
Supplier Catalog Number	GTX89196-PEP
Alternative Catalog Number	GTX89196-PEP-100
Manufacturer	GeneTex
Category	Proteine/Peptide
Species Reactivity	Mouse
Sensitivity	The peptide was used in the production of GTX89196.
NCBI	18706
Buffer	Reconstitute with 200 $\mu$ l distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H2O, no preservative.
Form	Lyophilized powder
Sequence	C-HQWLKDKNKGEIYD