

Diagnostica Vertrieb GmbH, Leipziger Straße 4 85386 Eching, Germany Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99 E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

GALR3 blocking peptide GTX89373-PEP

| Article Name | GALR3 blocking peptide |
|----------------------------|---|
| Biozol Catalog Number | GTX89373-PEP |
| Supplier Catalog Number | GTX89373-PEP |
| Alternative Catalog Number | GTX89373-PEP-100 |
| Manufacturer | GeneTex |
| Category | Proteine/Peptide |
| Species Reactivity | Human |
| Sensitivity | The peptide was used in the production of GTX89373. |
| NCBI | 8484 |
| Buffer | Reconstitute with 200 μ l distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H2O, no preservative. |
| Form | Lyophilized powder |
| Sequence | C-QGPEPREGPVHGGE |