



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## DARC blocking peptide GTX89467-PEP

Article Name	DARC blocking peptide
Biozol Catalog Number	GTX89467-PEP
Supplier Catalog Number	GTX89467-PEP
Alternative Catalog Number	GTX89467-PEP-100
Manufacturer	GeneTex
Category	Proteine/Peptide
Species Reactivity	Human
Sensitivity	The peptide was used in the production of GTX89467.
NCBI	2532
Buffer	Reconstitute with 200 $\mu$ l distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H2O, no preservative.
Form	Lyophilized powder
Sequence	HRAELSPSTENSSQLDFED-C
Application Notes	The antibodies(GTX89467) come at 0.5mg/ml and peptide as a 100ug pellet. When peptides are reconstituted in 200ul water, the concentration would be also 0.5mg/ml. To start, the best ratio would be 1:1 (which means molar excess of peptides relative to ant