



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Fluorescent Dye 488-I Streptavidin ABN-U0280

Article Name	Fluorescent Dye 488-I Streptavidin
Biozol Catalog Number	ABN-U0280
Supplier Catalog Number	U0280
Alternative Catalog Number	ABN-U0280-200
Manufacturer	Abnova
Category	Sonstiges
Application	Conjugation
Product Description	Fluorescence Dye 488-I with Streptavidin conjugate. Fluorescence Dye 488-I conjugates are significantly brighter than fluorescein conjugates and much more photostable. Additionally, the fluorescence of Fluorescence Dye 488-I is not affected by pH bet
Buffer	Water
Application Dilute	Labeling via Streptavidin for fluorescence imaging and other fluorescence-based biochemical analysis. For most fluorescent streptavidin applications, a concentration of 1-20 ug/mL is sufficient. The optimal working dilution should be determined by the