



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Fluorescent Dye 555-I Streptavidin ABN-U0288

Article Name	Fluorescent Dye 555-I Streptavidin
Biozol Catalog Number	ABN-U0288
Supplier Catalog Number	U0288
Alternative Catalog Number	ABN-U0288-200
Manufacturer	Abnova
Category	Sonstiges
Application	Conjugation
Product Description	Fluorescence Dye 555-I with Streptavidin conjugate. Compared with Cy3 conjugates, those of Fluorescence Dye 555-I are slightly redshifted, resulting in an optimal match to filters designed for Cy3 dyes in accompany with improved photostability. Besi
Buffer	Water
Application Dilute	Labeling via Streptavidin for fluorescence imaging and other fluorescence-based biochemical analysis. For most fluorescent streptavidin applications, a concentration of 1-20 ug/mL is sufficient. The optimal working dilution should be determined by the