



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Fluorescent Dye 594-I Streptavidin ABN-U0291

Article Name	Fluorescent Dye 594-I Streptavidin
Biozol Catalog Number	ABN-U0291
Supplier Catalog Number	U0291
Alternative Catalog Number	ABN-U0291-200
Manufacturer	Abnova
Category	Sonstiges
Application	Conjugation
Product Description	Fluorescence Dye 594-I with Streptavidin conjugate. Biomolecules conjugated to Fluorescence Dye 594-I exhibit little spectral overlap with green-fluorescent conjugates, and can be efficiently excited by 568nm line of Ar-Kr laser and by the 594nm line
Buffer	Water
Application Dilute	Labeling via Streptavidin for fluorescence imaging and other fluorescence-based biochemical analysis. For most fluorescent streptavidin applications, a concentration of 1-20 ug/mL is sufficient. The optimal working dilution should be determined by the