



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Fluorescent Dye 633-I Streptavidin ABN-U0295

Article Name	Fluorescent Dye 633-I Streptavidin
Biozol Catalog Number	ABN-U0295
Supplier Catalog Number	U0295
Alternative Catalog Number	ABN-U0295-200
Manufacturer	Abnova
Category	Sonstiges
Application	Conjugation, FC
Product Description	Fluorescence Dye 633-I with Streptavidin conjugate. Fluorescence emission of Fluorescence Dye 633-I is well separated from that of other commonly used red fluorophores, such as TAMRA, Texas Red, Alexa Fluor 594, Fluorescent Dye 594-I and R-phycoeryth
Buffer	Water
Application Dilute	Labeling via Streptavidin for fluorescence imaging and other fluorescence-based biochemical analysis. For most fluorescent streptavidin applications, a concentration of 1-20 ug/mL is sufficient. The optimal working dilution should be determined by the