



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Fluorescent Dye 700-I Streptavidin ABN-U0301

Article Name	Fluorescent Dye 700-I Streptavidin
Biozol Catalog Number	ABN-U0301
Supplier Catalog Number	U0301
Alternative Catalog Number	ABN-U0301-200
Manufacturer	Abnova
Category	Sonstiges
Application	Conjugation, FC
Product Description	Fluorescence Dye 700-I with Streptavidin conjugate. In comparison with Alexa Fluor 700, Fluorescence Dye 700-I are brighter with stronger absorption at 633nm. Fluorescence emission of Fluorescence Dye 700-I is well separated from that of other common
Buffer	Water
Application Dilute	Labeling via Streptavidin for fluorescence imaging and other fluorescence-based biochemical analysis. For most fluorescent streptavidin applications, a concentration of 1-20 ug/mL is sufficient. The optimal working dilution should be determined by the